

DHS/DCS Plate (Long/Short Barrel)

The Dynamic Hip Screw is designed to provide strong and stable internal fixation of a variety of intertrochanteric, subtrochanteric and basilar neck fractures, with minimal soft tissue irritation. The Dynamic Condylar Screw is designed to provide strong and stable internal fixation of certain distal femoral and subtrochanteric fractures, with minimal soft tissue irritation.




1. Preoperative planning

The size and angle of the plate as well as the length of the DHS Blade or Screw can be determined preoperatively by using the DHS Goniometer.

Note: If the DHS Blade is from 65 to 75 mm, a DHS plate with short barrel should be used to allow for Sufficient dynamization.

2. Position patient

Place the patient in a supine position on the operating table.



3. Reduce fracture

If possible, reduce the fracture under the image intensifier by means of traction, abduction and internal Rotation.

4. Access

Make a straight lateral skin incision of approximately 15 cm in length, starting two finger-widths proximal to the tip of the greater trochanter. Split the iliotibial tract lengthwise. Detach the m. vastus lateralis dorsally to the intramuscular membrane, retract ventrally and, if necessary, make a slight notch in the muscle in the -region of the innominate tubercle. Expose the proximal femoral shaft without retracting the periosteum.



DHS/DCS Plate Dynamic Hip and Condylar Screw

Indications and Contraindications

DHS

Including all combinations of DHS Screw, DHS Blade, DHS plate with DCP holes, LCP DHS plate and LCP DHS With collar-.


Indications DHS

- Pertrochanteric fractures of type 31-A1 and 31-A2
- Intertrochanteric fractures of type 31-A3
- Basilar neck fractures 31-B
(DHS Screw in conjunction with an antirotation screw)
- Subtrochanteric fractures

Contraindications DHS

The DHS is not to be used in cases where there is a high incidence of:

- Sepsis
 - Malignant primary or metastatic tumors
 - Material sensitivity
 - Compromised vascularity
- Recommendations DHS
- DHS Blade: for osteoporotic patients
 - DHS Screw B 14 mm: for revisions of DHS Screws B 13 mm
 - LCP DHS: for the use of shorter plates, especially in the case of femoral neck fractures
 - For certain subtrochanteric fractures, a 95° DCS plate is recommended.



31-B1

31-B2

31-B3


Trochanter Stabilizing Plate Indications

- Unstable pertrochanteric fractures of type 31-A2 and 31-A3, especially multi fragmentary fractures with a separated or longitudinally split greater trochanter

DCS

Indications DCS


- Proximal femur: Very proximally located, purely subtrochanteric fractures of types 32-A and 32-B
- Distal femur: Fractures of type 33-A (extra-articular, supracondylar) and fractures of type 33-C (fully articular Fractures)



33-A1

33-A2

33-A3



Contraindications DCS

- Pertrochanteric fractures or trochanteric fractures with Subtrochanteric expansion (31-A3)

33-C1

33-C2

33-C3

The Dynamic Hip Screw is designed to provide strong and stable internal fixation of a variety of intertrochanteric, subtrochanteric and basilar neck fractures, with minimal soft tissue irritation.

Strong

The DHS plates are made of 316L stainless steel and are cold worked strength.

Stable

The number of screw holes per plate length is maximized, without compromising plate's strength. This allows an increased number of fixation points with a smaller incision.

DCP (dynamic compression plate) holes in the DHS side plate:

- Allow angulations of 4.5 mm cortex screws, for lag screw fixation of medial fragments, and
- Allow axial compression and multiple-screw fixation of the main fragment in subtrochanteric fractures with shaft extension

Two flats within the DHS plate barrel correspond to the two-flat design of the DHS / DCS lag screw, preventing rotation of the lag screw within the barrel. The two-flat design also eases insertion of the plate over the DHS/DCS lag screw.



Minimal soft tissue irritation

The DHS plates have a low-profile design, reducing the risk of trochanteric bursitis.

The DHS plates are available in a wide range of sizes and barrel angles, with standard or short barrels, for varied clinical situations. The DHS/DCS lag screw, available from 50 mm to 145 mm lengths, easily glides within the DHS plate barrel for controlled collapse and impaction of fragments. When the fracture requires additional intraoperative compression, the DHS/DCS compression screw can be used; only one size compression screw is needed.

The DHS Instruments provide direct measurements throughout the DHS procedure, allowing proper reaming, tapping and lag screw insertion depth. The built-in stop and locking nut on the DHS triple reamer prevent over-reaming.

Dynamic Condylar Screw (DCS)

The Dynamic Condylar Screw is designed to provide strong and stable internal fixation of certain distal femoral and subtrochanteric fractures, with minimal soft tissue irritation.

Strong

The DCS plates are made of 316L stainless steel and are cold-worked for strength.

Stable

The two holes closest to the barrel accept 6.5 mm cancellous bone screws. This enhances stability by allowing additional fixation of the most distal condylar fracture fragments or fixation of the most proximal subtrochanteric fracture fragments.

DHS/DCS Plate Dynamic Hip and Condylar Screw

DCP holes in the DCS side plate allow angulation of 4.5 mm cortex screws and axial compression across a Shaft fracture.

The number of screw holes per plate length is maximized, without compromising plate strength. This allows an increased number of fixation points with a smaller incision.


Two flats within the DCS plate barrel correspond to the two-flat design of the lag screw, preventing rotation of the DHS/DCS lag screw within the barrel. The two-flat design also eases insertion of the plate over the DHS / DCS lag screw.

Minimal soft tissue irritation

The low-profile design reduces the risk of iliotibia band irritation (distal femoral fractures) and trochanteric bursitis (subtrochanteric fractures).

The DCS plates are available with 6 to 16 holes, for varied clinical situations. The DHS/DCS lag screw is available in 50 mm to 145 mm lengths. The DHS/DCS compression screw can be used for additional compression; only one size compression screw is needed.

The DCS instruments also provide direct measurements throughout the DCS procedure, allowing proper reaming, tapping, and lag screw insertion depth. The built-in stop and locking nut on the DCS triple reamer prevent over-reaming



Indications

DHS

The DHS is indicated for the following fractures of the proximal femur:

- Intertrochanteric fractures
- Subtrochanteric fractures*
- Basilar neck fractures

The DHS is indicated for stable fractures, and unstable fractures in which a stable medial buttress can be reconstructed. The DHS provides controlled collapse and compression of fracture fragments. This results in stable fixation and prevents undue stress concentration on the implant.

- For certain subtrochanteric fractures, a 95° device is the implant of choice (See “Using the DCS for Subtrochanteric Fractures,”).

DCS

The DCS is indicated for the following fractures of the distal femur:

- Intercondylar fractures**
- Supracondylar fractures**
- Unicondylar fractures**

** The following anatomic conditions should exist:

- 4 cm of distal femur should remain intact to provide support for the implant.
- A distal portion of the medial condyle should be intact for the DHS/DCS lag screw to gain good purchase.
If these conditions do not exist, a 4.5 mm LCP condylar plate or LCP distal femur plate should be considered.

The DCS is indicated for the following fractures of the proximal femur:

- Transverse subtrochanteric fractures[†]
- Short oblique subtrochanteric fractures[†]
- Long oblique subtrochanteric fractures[†]

[†] With the lesser trochanter avulsed or on the distal fragment (femoral shaft).¹

Plate selection

Barrel length

The standard 38 mm barrel length is most commonly indicated.

The 25 mm short barrel is indicated for specific clinical situations, including:

- Cases in which the standard barrel may not provide sufficient glide for the lag screw; i.e., a long impaction distance is expected
- A medial displacement osteotomy
- Unusually small femurs

Barrel angle

An evaluation of the angle subtended between the femoral neck and shaft axes (CCD, or collum-center-diaphysis angle) of the uninjured femur will aid in the selection of the most appropriate barrel angle. The 135° barrel angle is most commonly indicated.

Note: Greater barrel angles may produce biomechanical advantages in unstable cases; i.e., better gliding characteristics and reduced bending stresses on the plate/barrel junction, although correct placement of the implant becomes technically more difficult as barrel angles increase.²



SURGICAL TECHNIQUE

DHS PLATE

Insert a new DHS/DCS guide wire at the desired angle with the correct angled guide. The guide wire should be placed in the middle of the femoral head and should -extend into the subchondral bone.

Read the length of the DHS Screw or Blade directly off the guide wire with the measuring device.

If the guide wire is inserted into the subchondral bone -remove 10 mm from the measurement Assemble the triple reamer.

Slide the reamer over the drill bit until it clicks into place.

Set the triple reamer at the length of the implant Selected.

Mount the centering sleeve from the side onto the tap and lock it into place by turning the inner sleeve clock-wise against the outer sleeve.

Insert the connecting screw into the wrench, slide an appropriate DHS plate onto it and connect the DHS Screw to the wrench. For DHS screws shorter than or equal to 75 mm, take a DHS plate with short barrel. Mount the centering sleeve onto the wrench.

Insert the DHS guide shaft for coupling screw shaft, Insert the DHS impactor, Insert the DHS Compression Screw with Screw

Driver.



- 5400-INS-0014** D.H.S ANGLE GUIDE-135 °
- 5400-INS-0013** Ø 2.5 MM THREADED GUIDE WIRE -230 MM
- 5400-INS-0012** DIRECT MEASURING DEVICE ALUMINIUM
- 5400-INS-0002** D.H.S. TRIPLE REAMER Q. C.
- 5400-INS-0008** Ø 12.5 MM DHS TAP
- 5400-INS-0004** D.H.S/D.C.S CENTERING SLEEVE, LONG FOR T-WRENCH



- 5400-INS-0010** DHS/DCS GUIDE SHAFT FOR COUPLING SCREW SHORT
- 5400-INS-0009** DHS/DCS COUPLING SCREW, SHORT FOR INSERTING DHS SCREW
- 5400-INS-0011** D.H.S IMPACTOR
- 5400-INS-0016** Ø 4.5 MM SCREW DRIVER-9"

Use the drill guide and the drill bit to drill holes in a neutral position through the plate holes. Insert self tapping 4.5 mm cortex screws of appropriate length. Insert the 4.5 mm drill bit Q.C.



- 5400-INS-0026** Ø 4.5 MM/3.2 MM DOUBLE DRILL SLEEVE NORMAL
- 5400-INS-0019** Ø 4.5 MM DRILL BIT Q. C. (8")
- 5400-INS-0022** Ø 4.5 MM DEPTH GAUGE
- 5400-INS-0016** Ø 4.5 MM SCREW DRIVER-9"
- 5400-INS-0023** Ø 4.5 MM SCREW DRIVER WITH HOLDING SLEEVE -9"

Instruments

5400-INS-0001 Q. C. T-HANDLE



5400-INS-0002 D.H.S. TRIPLE REAMER Q. C



5400-INS-0003 D.C.S. TRIPLES REAMER Q. C.



5400-INS-0004 D.H.S/D.C.S CENTERING SLEEVE, LONG
FOR T-WRENCH



5400-INS-0005 D.H.S/D.C.S T-WRENCH



5400-INS-0006 DHS/DCS COUPLING SCREW, LONG FOR
DHS SCREW REMOVAL



5400-INS-0007 D.H.S/D.C.S CENTERING SLEEVE, SHORT
FOR DHS SCREW TAP



5400-INS-0008 \varnothing 12.5 MM DHS TAP



5400-INS-0009 DHS/DCS COUPLING SCREW, SHORT
FOR INSERTING DHS SCREW



5400-INS-0010 DHS/DCS GUIDE SHAFT FOR COUPLING
SCREW SHORT



5400-INS-0011 D.H.S IMPACTOR



5400-INS-0012 DIRECT MEASURING DEVICE ALUMINIUM



5400-INS-0013 Ø 2.5 MM THREADED GUIDE WIRE
230 MM



5400-INS-0014 D.H.S ANGLE GUIDE 135 °



5400-INS-0015 D.C.S. ANGLE GUIDE 95°



5400-INS-0016 Ø 4.5 MM SCREW DRIVER-9"



5400-INS-0017 Ø 4.5 MM BONE TAP



5400-INS-0018 Ø 3.2 MM DRILL BIT Q. C. (8")



5400-INS-0019 Ø 4.5 MM DRILL BIT Q. C. (8")



5400-INS-0020 Q. C. T-HANDLE FOR DRILL BIT USE



5400-INS-0021 VARIABLE ANGLE GUIDE



5400-INS-0022 Ø 4.5 MM DEPTH GAUGE



5400-INS-0023 4.5 MM SCREW DRIVER WITH HOLDING
SLEEVE -9"



5400-INS-0024 Ø 4.5 MM SCREW DRIVER QUICK COUPLING



5400-INS-0025 SCREW HOLDING FORCEP



5400-INS-0026 Ø 4.5 MM/3.2 MM DOUBLE DRILL
SLEEVE NORMAL



5400-INS-0027 Ø 4.5 MM DCP DRILL GUIDE
(NEUTRAL & LOADED)



5400-INS-0028 Ø 4.5 MM BONE TAP CORTICAL WITH T-HANDLE



5400-INS-0029 Ø 6.5 MM BONE TAP CANCELLOUS
WITH T-HANDLE



Implant Size:

125° LCP DHS PLATE LONG BARREL

SS316L	Titanium	Size
1369-SS-2503	1369-TT-2503	3 HOLES
1369-SS-2504	1369-TT-2504	4 HOLES
1369-SS-2505	1369-TT-2505	5 HOLES
1369-SS-2506	1369-TT-2506	6 HOLES
1369-SS-2507	1369-TT-2507	7 HOLES
1369-SS-2508	1369-TT-2508	8 HOLES
1369-SS-2509	1369-TT-2509	9 HOLES
1369-SS-2510	1369-TT-2510	10 HOLES
1369-SS-2512	1369-TT-2512	12 HOLES



130° LCP DHS PLATE LONG BARREL

SS316L	Titanium	Size
1369-SS-3003	1369-TT-3003	3 HOLES
1369-SS-3004	1369-TT-3004	4 HOLES
1369-SS-3005	1369-TT-3005	5 HOLES
1369-SS-3006	1369-TT-3006	6 HOLES
1369-SS-3007	1369-TT-3007	7 HOLES
1369-SS-3008	1369-TT-3008	8 HOLES
1369-SS-3009	1369-TT-3009	9 HOLES
1369-SS-3010	1369-TT-3010	10 HOLES
1369-SS-3012	1369-TT-3012	12 HOLES
1369-SS-3014	1369-TT-3014	14 HOLES
1369-SS-3016	1369-TT-3016	16 HOLES

135° LCP DHS PLATE LONG BARREL

SS316L	Titanium	Size
1369-SS-3503	1369-TT-3503	3 HOLES
1369-SS-3504	1369-TT-3504	4 HOLES
1369-SS-3505	1369-TT-3505	5 HOLES
1369-SS-3506	1369-TT-3506	6 HOLES
1369-SS-3507	1369-TT-3507	7 HOLES
1369-SS-3508	1369-TT-3508	8 HOLES
1369-SS-3509	1369-TT-3509	9 HOLES
1369-SS-3510	1369-TT-3510	10 HOLES
1369-SS-3512	1369-TT-3512	12 HOLES
1369-SS-3514	1369-TT-3514	14 HOLES
1369-SS-3516	1369-TT-3516	16 HOLES

140° LCP DHS PLATE LONG BARREL

SS316L	Titanium	Size
1369-SS-4003	1369-TT-4003	3 HOLES
1369-SS-4004	1369-TT-4004	4 HOLES
1369-SS-4005	1369-TT-4005	5 HOLES
1369-SS-4006	1369-TT-4006	6 HOLES
1369-SS-4007	1369-TT-4007	7 HOLES
1369-SS-4008	1369-TT-4008	8 HOLES
1369-SS-4009	1369-TT-4009	9 HOLES
1369-SS-4010	1369-TT-4010	10 HOLES
1369-SS-4012	1369-TT-4012	12 HOLES
1369-SS-4014	1369-TT-4014	14 HOLES
1369-SS-4016	1369-TT-4016	16 HOLES

125° LCP DHS PLATE SHORT BARREL

SS316L	Titanium	Size
1370-SS-2503	1370-TT-2503	3 HOLES
1370-SS-2504	1370-TT-2504	4 HOLES
1370-SS-2505	1370-TT-2505	5 HOLES
1370-SS-2506	1370-TT-2506	6 HOLES
1370-SS-2507	1370-TT-2507	7 HOLES
1370-SS-2508	1370-TT-2508	8 HOLES
1370-SS-2509	1370-TT-2509	9 HOLES
1370-SS-2510	1370-TT-2510	10 HOLES
1370-SS-2512	1370-TT-2512	12 HOLES



130° LCP DHS PLATE SHORT BARREL

SS316L	Titanium	Size
1370-SS-3003	1370-TT-3003	3 HOLES
1370-SS-3004	1370-TT-3004	4 HOLES
1370-SS-3005	1370-TT-3005	5 HOLES
1370-SS-3006	1370-TT-3006	6 HOLES
1370-SS-3007	1370-TT-3007	7 HOLES
1370-SS-3008	1370-TT-3008	8 HOLES
1370-SS-3009	1370-TT-3009	9 HOLES
1370-SS-3010	1370-TT-3010	10 HOLES
1370-SS-3012	1370-TT-3012	12 HOLES
1370-SS-3014	1370-TT-3014	14 HOLES
1370-SS-3016	1370-TT-3016	16 HOLES

135° LCP DHS PLATE SHORT BARREL

SS316L	Titanium	Size
1370-SS-3503	1370-TT-3503	3 HOLES
1370-SS-3504	1370-TT-3504	4 HOLES
1370-SS-3505	1370-TT-3505	5 HOLES
1370-SS-3506	1370-TT-3506	6 HOLES
1370-SS-3507	1370-TT-3507	7 HOLES
1370-SS-3508	1370-TT-3508	8 HOLES
1370-SS-3509	1370-TT-3509	9 HOLES
1370-SS-3510	1370-TT-3510	10 HOLES
1370-SS-3512	1370-TT-3512	12 HOLES
1370-SS-3514	1370-TT-3514	14 HOLES
1370-SS-3516	1370-TT-3516	16 HOLES

DHS/DCS Plate Dynamic Hip and Condylar Screw

140° LCP DHS PLATE SHORT BARREL

SS316L	Titanium	Size
1370-SS-4003	1370-TT-4003	3 HOLES
1370-SS-4004	1370-TT-4004	4 HOLES
1370-SS-4005	1370-TT-4005	5 HOLES
1370-SS-4006	1370-TT-4006	6 HOLES
1370-SS-4007	1370-TT-4007	7 HOLES
1370-SS-4008	1370-TT-4008	8 HOLES
1370-SS-4009	1370-TT-4009	9 HOLES
1370-SS-4010	1370-TT-4010	10 HOLES
1370-SS-4012	1370-TT-4012	12 HOLES
1370-SS-4014	1370-TT-4014	14 HOLES
1370-SS-4016	1370-TT-4016	16 HOLES

95° LCP DCS PLATE

SS316L	Titanium	Size
1322-SS-9504	1322-TT-9504	4 HOLES
1322-SS-9505	1322-TT-9505	5 HOLES
1322-SS-9506	1322-TT-9506	6 HOLES
1322-SS-9507	1322-TT-9507	7 HOLES
1322-SS-9508	1322-TT-9508	8 HOLES
1322-SS-9509	1322-TT-9509	9 HOLES
1322-SS-9510	1322-TT-9510	10 HOLES
1322-SS-9511	1322-TT-9511	11 HOLES
1322-SS-9512	1322-TT-9512	12 HOLES
1322-SS-9514	1322-TT-9514	14 HOLES
1322-SS-9516	1322-TT-9516	16 HOLES
1322-SS-9518	1322-TT-9518	18 HOLES



Screws:

Ø4.5 MM CORTICAL SCREW - 14TPI - SELF TAPPING

1111-SS-4510/1111-SS-4580 SS 10 mm to 80 mm

1111-TT-4510/1111-TT-4580 TT 10 mm to 80 mm



Ø12.5 MM DHS SCREW WITH COMPRESSION SCREW

1123-SS-1245/ 1123-SS-12120 SS 45 mm to 120 mm

1123-TT-1245/ 1123-TT-12120 TT 45 mm to 120 mm



Ø6.5 MM LCP CANCELLOUS SCREW 16 MM THREADED

1421-SS-6525/1421-SS-65120 SS 25 mm to 120 mm

1421-TT-6525/1421-TT-65120 TT 25 mm to 120 mm



Ø6.5 MM LCP CANCELLOUS SCREW 32 MM THREADED

1422-SS-6540/1422-SS-65120 SS 40 mm to 120 mm

1422-TT-6540/1422-TT-65120 TT 40 mm to 120 mm



Ø6.5 MM LCP CANCELLOUS SCREW FULLY THREADED

1423-SS-6525/1423-SS-65100 SS 25 mm to 100 mm

1423-TT-6525/1423-TT-65100 TT 25 mm to 100 mm



Address:

Corporate Office

310, Sanket Avenue, Opp. Ambe Vidhyalaya,
Sama-savli road Vadodara-390024

Plant Office:

R.S.No. 99/2, P-2, Village & Post Karmasiya
Kheda Rasulabad – Vitoj Road,
Taluka - Waghodiya Dist. Vadodara-391510
Gujarat –India.